



Brewer
more possibilities

ACCESS HIGH-LOW EXAM TABLE

BrewerAccess High-Low


EXAM TABLE Designed to take patients to a new height of care.

The Brewer Access High-Low

Exam Table is clearly the:

- most convenient for patient consultations and examinations
- safest for patient and wheelchair transfers
- most reliable
- most user-friendly
- most patient-friendly
- most technologically advanced high-low exam table on the market

**See why Brewer means...
more possibilities.**



Chair



Table



Chair-to-table design effortlessly takes patient from consultation to examination
18" low height for easy and safe patient transfers
Return to chair feature easily takes table back to chair position in one, simple step (on power back model only)



Standard grab bars provide patient added safety and comfort
Provides patients a stable and secure grab point to use when getting on and off table
Conveniently swings out for lateral wheelchair transfers
Easily removable
Standard on all models

Abundant storage capacity
Patented pass-through side drawers allow left or right access
Eliminates right or left hand setup problems
Additional spacious front drawer
Convenient covered paper roll storage holds up to six rolls
Cover keeps paper rolls clean



Pass-through work surface pulls out on right or left side
Saves space by eliminating the need for additional stands or trays
Perfect for charts, instruments, procedure preparation, or portable computer



Comfortable, ergonomic upholstery
Easy to clean surface; seamless design



Intuitive foot control
On / off switch (located on the back of table) allows table's foot control to be disabled and prevents patients and / or children from moving table



Extendable, adjustable ergonomic stirrups

Brewer's unique approach to product development is centered on researching the needs of healthcare professionals and their patients.

Before a design begins, we spend countless hours talking to doctors, nurses, and patients to understand what is most important to them.

The result is a superior product that delivers the exact features needed to transform the exam room into a safe, efficient, and comfortable environment. We invite you to discover the difference...

BrewerAccess High-Low

EXAM TABLE Designed to take patients to a new height of care.

Durability

- 3-year warranty

Safety

- 18" low height
- Easy on / off for patients
- Safety grab bars for patient stability and positioning
- 450 lb. weight capacity*
- Safety, extra low voltage
- Standard on / off switch

Function

- Patented pass-through side drawers
- Largest storage capacity available
- Pass-through work surface on side of table
- Backrest adjusts for easy patient positioning (pneumatic or power back available)
- Additional paper storage
- Return to chair feature (on power back model)
- Easy to clean, seamless upholstery
- Standard foot control

Patient Comfort

- Easy lateral wheelchair transfers
- Safety grab bars
- Multi-density support upholstery
- Adjustable ergonomic stirrups
- Toe board / leg extension / plastic debris tray
- Premium, plush upholstery available

Aesthetics

- 18 standard upholstery color options
- 6 premium, plush upholstery color options
- Custom colors available

Weight Capacity

- 450 lbs.*

Table Dimensions

- Height: 18" - 37"
- Length: 57" extends to 72"
- Upholstery width: 28"
- Leg extension length: 15"

Covered Paper Roll Storage / Paper Roll Holder

- Storage capacity: 23" x 23" x 3.25"
- Storage for up to 6 additional paper rolls
- Accepts 21" x 3" and 18" x 3" paper rolls

Storage

- Front drawer: 19.7" x 19.2" x 4"
- Side drawers: 22.7" x 20.2" x 7"
- Pass-through side drawers accessible from right or left side of table
- Covered paper roll storage

Certifications

- UL listed, CUL listed

Options

- Power back / return to chair
- Premium plush upholstery
- 5° pelvic tilt / 60-Watt drawer heater in front top drawer

#6000	pneumatic back
#6001	pneumatic back with pelvic tilt and drawer heater
#6500	power back, return to chair
#6501	power back, return to chair with pelvic tilt and drawer heater

accessories

#100407	front drawer dividers
#98464	stainless steel treatment pan
#100220	urology drain pan system
#99504	deluxe articulating knee crutches
#100230	deluxe articulating knee crutch clamp

*under normal use

All Brewer Access High-Low Exam Tables come with a 3-year warranty.

The Brewer Company reserves the right to improve or modify the design of its products consistent with new methods and techniques.

Medical Device Depot
Inc.

www.medicaldevicedepot.com
1-877-646-3300

Brewer
more possibilities