

# Lower Extremity Physiologic Study Multiple Level (Segmental Exam)

Patient Name: SEG PVR EXAM

Patient ID: 1233356

Age: 67

Exam Date: 06/13/08 12:54:39

Risk Factors	
<input checked="" type="checkbox"/> Tobacco Use	<input type="checkbox"/> Hypertension
<input checked="" type="checkbox"/> Diabetes	<input type="checkbox"/> Hyperlipidemia
<input type="checkbox"/> Heart Disease	<input type="checkbox"/> Stroke/TIA
<input type="checkbox"/> Previous CV Event	<input type="checkbox"/> Prev. Vascular Surgery
3 Pack per Day Smoker	


Current Symptoms	
<input type="checkbox"/> Intermittent Claudication	<input type="checkbox"/> Location
<input checked="" type="checkbox"/> Numbness, Tingling in Feet	<input type="checkbox"/> Buttocks
<input type="checkbox"/> Ulcerations	<input type="checkbox"/> Thigh
<input checked="" type="checkbox"/> Rest Pain	<input type="checkbox"/> Calf
<input type="checkbox"/> Gangrene	<input checked="" type="checkbox"/> Feet
No other symptoms	

ABI Severity Index for Ankles Only	
1.00 - 1.29	Normal
0.91 - 0.99	Borderline
0.41 - 0.90	Moderate/Mild
0.00 - 0.40	Severe


Comments: Exam performed by Sheila Smith

## RIGHT


Brachial  
120




Thigh  
123  
mmHg  
1.02  
Index




Abv Knee  
126  
mmHg  
1.05  
Index




Calf  
129  
mmHg  
1.07  
Index



PT  
132  
mmHg  
1.10  
Index




Toe  
98  
mmHg  
0.81  
Index




## LEFT


Brachial  
118




Thigh  
125  
mmHg  
1.04  
Index




Abv Knee  
126  
mmHg  
1.05  
Index



Calf  
129  
mmHg  
1.07  
Index



PT  
132  
mmHg  
1.10  
Index



Toe  
105  
mmHg  
0.87  
Index

Interpretation: Patient has normal waveform and pressure  
no PAD